

# Monitoring Fatty Liver Disease During Pre/Post-Bariatric Surgery With Non-Invasive LIVERFAST™

**Teresa Gonzalo, PharmD MBA, Medical Science Director |** Fibronostics, Inc.

**Mona Munteanu, MD, Senior Medical Affairs |** Fibronostics, Inc.


**Ronald Quiambao, MD, Chief Medical Officer |** Fibronostics, Inc.

**Roni Amiel, Chief Executive Officer |** Fibronostics, Inc.

## LIVER DISEASE AND BARIATRIC SURGERY

According to the American Society for Metabolic and Bariatric Surgery (ASMBS)<sup>1</sup>, non-alcoholic fatty liver disease (NAFLD) is one of the obesity-related co-morbidities that qualifies a patient to undergo a bariatric surgery if BMI > 35. By 2030, it is predicted that nearly half of adults in the USA will have obesity<sup>2</sup>. Over 80% of the patients with obesity submitted to bariatric surgery suffer from (NAFLD), with 25% - 55% resulting in steatohepatitis (NASH) and 2% - 12% liver fibrosis and cirrhosis<sup>3</sup>. Current management of obese patients with NASH consists of lifestyle recommendations as very few therapeutic strategies for NASH are available<sup>4</sup>. Unfortunately, lifestyle interventions, rarely permit more than 10% total body weight loss, the threshold for meaningful improvement in patients with NASH to reduce inflammation and fibrosis<sup>5</sup>.

Numerous clinical practice guidelines including AASLD, EASL-EASD-EASO, APASL, and WHO recommend non-invasive serum biomarker-based diagnostic modalities to diagnose NAFLD<sup>6-8</sup>.



## CPT Code: 0166U

**THE ONLY NON-INVASIVE BIOMARKER TECHNOLOGY APPROVED BY: AMA & CMS**

<b>Complete Staging of the Liver</b> Fibrosis Inflammation Steatosis	<b>Stage Liver injury for Common Liver Diseases</b> HBV HCV NAFLD ALD
<b>Non-Invasive &amp; Cost Effective</b> Easily repeatable for frequent follow-up & monitoring	<b>10 Biomarkers including</b> Alpha-2-Macroglobulin Haptoglobin Apolipoprotein A1

## LIVERFAST™: AN ADVANCED BLOOD BASED TEST FOR LIVER DISEASES

---

Using the latest available artificial intelligence (AI) technology based on neuronal networking, an advanced algorithm LIVERFAST™ has been developed on multiethnic population as non-invasive clinical tool for staging and grading NAFLD<sup>9</sup>. LIVERFAST™ individual serum biomarkers have been identified as appropriate biomarkers for liver disease evaluation and results from FDA cleared assays<sup>10</sup>.

LIVERFAST™ combines anthropometric (age, gender and BMI) and validated serum biomarkers including alpha-2-macroglobulin, haptoglobin, apolipoprotein A1, bilirubin, gamma glutamyl transpeptidase (GGT), alanine aminotransferase (ALT), aspartate aminotransferase (AST), triglycerides, cholesterol (total) and fasting glucose to generate a score and a categorization of steatohepatitis using similar to histology grading and staging of the elementary liver lesions<sup>11-15</sup>.

Clinicians overly rely on abnormal liver transaminases alone to identify patients with NAFLD and NASH, although they have insufficient reliability to detect fibrosis<sup>16,17</sup>.

LIVERFAST™ Fibrosis score has been greatly optimized that discriminate fibrosis from other liver lesions that occur concomitantly as steatosis and inflammation which is biased by fluctuations in ALT.

LIVERFAST™ Activity score detects the degree of ballooning and lobular inflammation and discriminates from the fibrosis and steatosis liver lesions.

LIVERFAST™ Steatosis score detects lower amounts of steatosis than liver ultrasound (5%) permitting the diagnosis of NAFLD and the monitoring of patients.

To reliably diagnose NASH, a simplified scoring system, the SAF score (steatosis, activity, fibrosis) intended for pathologists has been validated in patients with morbid obesity undergoing bariatric surgery<sup>18,19</sup>, LIVERFAST™ is using the cutoffs tightly adapted to NAFLD and as discriminative as the histological SAF classification for steatosis, activity and fibrosis for the non-invasive determination of NAFLD/NASH with a cost-efficient approach.

## NAFLD AND NASH DIAGNOSIS PRE AND POST-BARIATRIC SURGERY

---

The subsequent risk of cirrhosis and hepatocellular carcinoma (HCC) emphasizes an urgent need for effective therapy to reverse fibrosis in morbidly obese patients diagnosed with NASH<sup>20</sup>.

Bariatric surgery efficacy has been proven for achieving sustained weight loss in NAFLD and NASH patients and can reverse risk factors that contribute to the pathogenesis of NAFLD, including dyslipidemia, insulin resistance, and hepatic inflammation, making it a promising treatment option for NAFLD<sup>21-23</sup>.

Various studies evaluating the histologic impact of bariatric surgery reported the disappearance of NASH in over 80% of cases<sup>24-26</sup>.

However, the initial severity of liver disease may affect bariatric surgery outcomes<sup>27</sup>, as some

studies have shown that NASH-diagnosed patients have a higher risk of liver-related mortality in the long-term follow-up<sup>28</sup>. Therefore, screening liver damage may help in the decision-making process of high-risk patients.

### **PRE-BARIATRIC SURGERY: STRATIFYING PATIENT RISK WITH LIVERFAST™**

---

To prevent unwanted post-bariatric surgery events, it is crucial to individualize patient selection and management.

A challenging element of the diagnostic workup of patients with morbid obesity is the NAFLD determination of disease severity. The goals here are double: to identify NASH, among the factors proceeding with bariatric surgery, and to identify those with advanced NAFLD at increased morbidity and mortality.

LIVERFAST™ can evaluate liver damage by stratifying high-risk patients with well-known metabolic comorbidities associated with NAFLD like obesity, type 2 diabetes, hyperlipidemia, hypertension, and metabolic syndrome<sup>4</sup>. LIVERFAST™ simplifies the indication of bariatric

surgery by identifying liver diseases and NASH, the final report with patient's stage of liver fibrosis, activity and steatosis will be revised by the physician and communicated appropriately to the patient.

### **POST-BARIATRIC SURGERY: MONITORING PATIENTS WITH LIVERFAST™**

---


Despite initial concerns about rapid weight loss from bariatric surgery could exacerbate NASH or acute liver failure in morbidly obese patients<sup>29,30</sup>, more recent surgical techniques such as Roux-en-Y gastric bypass<sup>31,32</sup> and intragastric balloon<sup>33</sup> have shown improvements in liver histopathologic scoring after 5 years of follow-up evaluation.

Liver biopsy is not adapted to monitoring of high amounts of patients undergoing bariatric surgery.

Besides sample variability, liver biopsy has non-negligible morbidity and mortality. Therefore, LIVERFAST™ provides a surrogate of liver biopsy without risks and bypassing its limits.

LIVERFAST™, a quantitative score, provides more granularity than the categorical liver biopsy to evaluate post-bariatric surgery regression of NASH or even worsening that could be noted in a small subset of patients<sup>34-37</sup>.

LIVERFAST™ simplifies liver monitoring in post-bariatric patients being crucial for the appropriate management and prevention of worse outcomes.

**NAFLD and NASH Diagnosis Pre and Post-Bariatric Surgery**


The LIVERFAST™ report provides scores along with staging of liver disease. The results will guide the physician's medical decisions and, communicated appropriately to the patient, will improve its compliance.

Furthermore, LIVERFAST is the only test able to discriminate fibrosis from steatosis and activity without bias in fibrosis estimation related to the presence of activity or steatosis.

Non-invasive diagnostic tools such as LIVERFAST™ are easy to perform, less expensive, and readily available and aid to the early diagnosis and better prognosis in patients with NAFLD and NASH.

**ORDERING LIVERFAST™ FOR PATIENTS**

A simple blood draw will be performed to obtain the LIVERFAST™ using Fibronostics technology.

In order to obtain a LIVERFAST™ digital assay test report, a healthcare provider must prescribe the LIVERFAST™ and order it through the laboratory with the CPT code. LIVERFAST™ ordering is integrated with laboratory ordering systems and communicate the results to the physician. The serum biomarkers are obtained from a CLIA certified lab and the results are submitted to the LIVERFAST™ algorithm platform.

**Staging/Grading of three lesions – Fibrosis Activity & Steatosis**

	LIVERFAST™	Transient Elastography	MRI	Ultrasound	APRI	ELF	FIB-4
<b>Fibrosis</b>	YES	YES	YES	NO	YES	YES	YES
<b>Activity</b>	YES	NO	NO	NO	NO	NO	NO
<b>Steatosis</b>	YES	YES	YES	YES	NO	NO	NO
<b># Biomarkers</b>	10	NA	NA	NA	2	3	4

**LIVERFAST™ provides complete evaluation and staging of fibrosis, activity and steatosis**

(MRI = Magnetic Resonance Imaging, APRI = AST-to-Platelet Ratio Index, ELF = Enhanced Liver Fibrosis, FIB-4 = Index for Liver fibrosis)

23.  
24.

## References

- Clinical Practice Guidelines for the Support of the Bariatric Surgery Patient | ASMB. American Society for Metabolic and Bariatric Surgery. Published March 1, 2013. Accessed June 22, 2020. <https://asmbs.org/resources/clinical-practice-guidelines-for-the-perioperative-nutritional-metabolic-and-nonsurgical-support-of-the-bariatric-surgery-patient>
- Ward ZJ, Bleich SN, Cradock AL, et al. Projected U.S. State-Level Prevalence of Adult Obesity and Severe Obesity. *New England Journal of Medicine*. 2019;381(25):2440-2450. doi:10.1056/NEJMsa1909301
- Morita S, Neto DDS, Morita FHA, Morita NK, Lobo SMA. Prevalence of Non-alcoholic Fatty Liver Disease and Steatohepatitis Risk Factors in Patients Undergoing Bariatric Surgery. *Obes Surg*. 2015;25(12):2335-2343. doi:10.1007/s11695-015-1696-5
- Chalasan N, Younossi Z, Lavine JE, et al. The diagnosis and management of non-alcoholic fatty liver disease: practice Guideline by the American Association for the Study of Liver Diseases, American College of Gastroenterology, and the American Gastroenterological Association. *Hepatology*. 2012;55(6):2005-2023. doi:10.1002/hep.25762
- Promrat K, Kleiner DE, Niemeier HM, et al. Randomized controlled trial testing the effects of weight loss on nonalcoholic steatohepatitis. *Hepatology*. 2010;51(1):121-129. doi:10.1002/hep.23276
- AASLD-IDSA HCV Guidance Panel. Hepatitis C Guidance 2018 Update: AASLD-IDSA Recommendations for Testing, Managing, and Treating Hepatitis C Virus Infection. *Clin Infect Dis*. 2018;67(10):1477-1492. doi:10.1093/cid/ciy585
- European Association for the Study of the Liver (EASL), European Association for the Study of Diabetes (EASD), European Association for the Study of Obesity (EASO). EASL-EASD-EASO Clinical Practice Guidelines for the management of non-alcoholic fatty liver disease. *J Hepatol*. 2016;64(6):1388-1402. doi:10.1016/j.jhep.2015.11.004
- WHO | Guidelines for the screening, care and treatment of persons with chronic hepatitis C infection. WHO. Accessed October 23, 2019. <http://www.who.int/hepatitis/publications/hepatitis-c-guidelines-2016/en/>
- Aravind A, Bahrivani AG, Quiambao R, Gonzalo T. Machine Learning Technology for Evaluation of Liver Fibrosis, Inflammation Activity and Steatosis (LIVERFAST™). *Journal of Intelligent Learning Systems and Applications*. 2020;12(2):31-49. doi:10.4236/jilsa.2020.122003
- Zhou J-H, Cai J-J, She Z-G, Li H-L. Noninvasive evaluation of nonalcoholic fatty liver disease: Current evidence and practice. *World J Gastroenterol*. 2019;25(11):1307-1326. doi:10.3748/wjg.v25.i11.1307
- Atanasova E, Martinova F, Jelev D, et al. Alpha-2 Macroglobulin Is The Simplest Serum Biomarker For Liver Fibrosis And Fibrogenesis In Chronic Hepatitis C. *Journal of Medical and Dental Practice*. 2015;2:153-164. doi:10.18044/MedInform.201522.153
- Morling JR, Guha IN. Biomarkers of liver fibrosis. *Clinical Liver Disease*. 2016;7(6):139-142. doi:10.1002/cld.555
- Neuman MG, Cohen LB, Nanau RM. Biomarkers in nonalcoholic fatty liver disease. *Can J Gastroenterol Hepatol*. 2014;28(11):607-618. Accessed October 24, 2019. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4277175/>
- Wong VW-S, Adams LA, de Lédinghen V, Wong GL-H, Sookoian S. Noninvasive biomarkers in NAFLD and NASH - current progress and future promise. *Nat Rev Gastroenterol Hepatol*. 2018;15(8):461-478. doi:10.1038/s41575-018-0014-9
- Fallatah HI. Noninvasive Biomarkers of Liver Fibrosis: An Overview. *Advances in Hepatology*. 2014;2014:1-15. doi:10.1155/2014/357287
- Mofrad P, Contos MJ, Haque M, et al. Clinical and histologic spectrum of nonalcoholic fatty liver disease associated with normal ALT values. *Hepatology*. 2003;37(6):1286-1292. doi:10.1053/jhep.2003.50229
- de Lédinghen V, Vergniol J, Gonzalez C, et al. Screening for liver fibrosis by using FibroScan® and FibroTest in patients with diabetes. *Dig Liver Dis*. 2012;44(5):413-418. doi:10.1016/j.dld.2011.12.005
- Bedossa P. Utility and appropriateness of the fatty liver inhibition of progression (FLIP) algorithm and steatosis, activity, and fibrosis (SAF) score in the evaluation of biopsies of nonalcoholic fatty liver disease. *Hepatology*. 2014;60(2):565-575. doi:10.1002/hep.27173
- Bedossa P, Poitou C, Veyrie N, et al. Histopathological algorithm and scoring system for evaluation of liver lesions in morbidly obese patients. *Hepatology*. 2012;56(5):1751-1759. doi:10.1002/hep.25889
- Wadden TA, Volger S, Sarwer DB, et al. A Two-Year Randomized Trial of Obesity Treatment in Primary Care Practice. *New England Journal of Medicine*. 2011;365(21):1969-1979. doi:10.1056/NEJMoa1109220
- Mattar SG, Velcu LM, Rabinovitz M, et al. Surgically-Induced Weight Loss Significantly Improves Nonalcoholic Fatty Liver Disease and the Metabolic Syndrome. *Ann Surg*. 2005;242(4):610-620. doi:10.1097/01.sla.0000179652.07502.3f
- Weiner RA. Surgical Treatment of Non-Alcoholic Steatohepatitis and Non-Alcoholic Fatty Liver Disease. *DDI*. 2010;28(1):274-279. doi:10.1159/000282102
- Laursen TL, Hagemann CA, Wei C, et al. Bariatric surgery in patients with non-alcoholic fatty liver disease - from pathophysiology to clinical effects. *World J Hepatol*. 2019;11(2):138-149. doi:10.4254/wjh.v11.i2.138
- Barker KB, Palekar NA, Bowers SP, Goldberg JE, Pulcini JP, Harrison SA. Non-alcoholic steatohepatitis: effect of Roux-en-Y gastric bypass surgery. *Am J Gastroenterol*. 2006;101(2):368-373. doi:10.1111/j.1572-0241.2006.00419.x
- Schneck A-S, Anty R, Patouraux S, et al. Roux-En Y Gastric Bypass Results in Long-Term Remission of Hepatocyte Apoptosis and Hepatic Histological Features of Non-alcoholic Steatohepatitis. *Front Physiol*. 2016;7:344. doi:10.3389/fphys.2016.00344
- Klein S, Mittendorfer B, Eagon JC, et al. Gastric bypass surgery improves metabolic and hepatic abnormalities associated with nonalcoholic fatty liver disease. *Gastroenterology*. 2006;130(6):1564-1572. doi:10.1053/j.gastro.2006.01.042
- Elnahas A, Nguyen GC, Okrainec A, Queresby F, Jackson TD. The effect of underlying liver disease on short-term outcomes following bariatric surgery. *Surg Endosc*. 2014;28(9):2708-2712. doi:10.1007/s00464-014-3532-8
- Goossens N, Hoshida Y, Song WM, et al. Nonalcoholic Steatohepatitis Is Associated With Increased Mortality in Obese Patients Undergoing Bariatric Surgery. *Clinical Gastroenterology and Hepatology*. 2016;14(11):1619-1628. doi:10.1016/j.cgh.2015.10.010
- Luyckx FH, Lefebvre PJ, Scheen AJ. Non-alcoholic steatohepatitis: association with obesity and insulin resistance, and influence of weight loss. *Diabetes Metab*. 2000;26(2):98-106.
- Vargas V, Allende H, Lecube A, et al. Surgically induced weight loss by gastric bypass improves non alcoholic fatty liver disease in morbid obese patients. *World J Hepatol*. 2012;4(12):382-388. doi:10.4254/wjh.v4.i12.382
- Caiazzo R, Lassailly G, Leteurtre E, et al. Roux-en-Y gastric bypass versus adjustable gastric banding to reduce nonalcoholic fatty liver disease: a 5-year controlled longitudinal study. *Ann Surg*. 2014;260(5):893-898; discussion 898-899. doi:10.1097/SLA.0000000000000945
- Mathurin P, Hollebecque A, Arnalsteen L, et al. Prospective study of the long-term effects of bariatric surgery on liver injury in patients without advanced disease. *Gastroenterology*. 2009;137(2):532-540. doi:10.1053/j.gastro.2009.04.052
- Bazerbachi F, Vargas EJ, Rizk M, et al. Intra-gastric Balloon Placement Induces Significant Metabolic and Histologic Improvement in Patients With Nonalcoholic Steatohepatitis. *Clin Gastroenterol Hepatol*. Published online April 30, 2020. doi:10.1016/j.cgh.2020.04.068
- Silverman EM, Sapala JA, Appelman HD. Regression of hepatic steatosis in morbidly obese persons after gastric bypass. *Am J Clin Pathol*. 1995;104(1):23-31. doi:10.1093/ajcp/104.1.23
- Moretto M, Kupski C, da Silva VD, Padoin AV, Mottin CC. Effect of bariatric surgery on liver fibrosis. *Obes Surg*. 2012;22(7):1044-1049. doi:10.1007/s11695-011-0559-y
- Schwenger KJP, Fischer SE, Jackson T, Okrainec A, Allard JP. In nonalcoholic fatty liver disease, Roux-en-Y gastric bypass improves liver histology while persistent disease is associated with lower improvements in waist circumference and glycemic control. *Surg Obes Relat Dis*. 2018;14(9):1233-1239. doi:10.1016/j.soard.2018.06.007
- Lee Y, Doumouras AG, Yu J, et al. Complete Resolution of Nonalcoholic Fatty Liver Disease After Bariatric Surgery: A Systematic Review and Meta-analysis. *Clin Gastroenterol Hepatol*. 2019;17(6):1040-1060.e11. doi:10.1016/j.cgh.2018.10.017